

DANI

UNSOPHISTICATION IS THE ULTIMATE SOPHISTICATION

Master
Master AS 2

AUTOMATED DUAL SIMULTANEOUS INJECTION

UTMOST PRODUCTIVITY AND SELECTIVITY
THROUGH SIMPLICITY

The image shows a Master AS 2 Automated Dual Injection autosampler. It is a tall, white, rectangular machine with a central vertical column. At the bottom, there is a sample tray containing a grid of small, blue-capped vials. The machine has a sleek, modern design with various ports and a control panel on the left side. The background is a light gray with decorative blue and white geometric shapes.

Master AS 2

Automated Dual Injection

Master AS 2 Automated DUAL SIMULTANEOUS Injection

Dual Simultaneous Injection into two capillary columns at the same time


Autosamplers can significantly augment **Result Repeatability** and **Extend Analytical Capabilities**.

The Master AS 2 Double Injection provides improved **Productivity** and **Flexibility**.

Master AS 2
Automated
DUAL SIMULTANEOUS Injection

Based on DANI Master AS, the Master AS 2 Dual Injection is a robotic X-Y-Z coordinate sampling system providing unique simultaneous injection into two equal or different capillary columns.

Increased sample throughput to meet each laboratory's changing workload.



DANI

The Master AS Liquid Autosampler and Master AS2 Double Injection allow to process vials sequentially and unattended with enhanced precision and accuracy, providing increased sample throughput and decreased cost per sample.


The easy-to-use liquid autosamplers are totally controlled through the Master GC touchscreen. Analysis can be performed unattended by using the straightforward interface to set all parameters and the AS sequence.

Unlike built-in rack autosamplers, the extremely flexible **Master AS** and **Master AS 2** autosamplers feature three separate and removable 2-mL vial racks for a total capacity of 160 samples.

10mL solvent and waste vials can be used allowing numerous solvent possibilities to eliminate carry-over. Additionally, the sample trays are delocalized from the GC oven, preventing possible sample degradation and solvent evaporation caused by heat exposure.

All injection parameters can be optimized according to the laboratory's needs, the sample type, and the application. Up to seven different syringe capacity from 5 μ L to 500 μ L can be used offering a great choice of injection volumes.

The Master AS autosamplers present unmatched injection capabilities. The same vial can be sampled up to 100 times and the sample can be introduced into three different injectors without requiring autosampler re-alignment or any calibration procedure.



Automation for Accurate Results

Flexibility & Productivity

Double Productivity : through two trays of 80 2mL vials each

Wide range of syringes available : 5, 10, 25, 50, 100, 250, 500 μ L

Injection volume range : 0.1 μ L - 500 μ L

Priority Sequence

Internal Standard Addition

The Master AS Double Injection provide increased Precision and Accuracy through several key implementations :

Parameter Controls:

Pre- and post- injection solvent washing
Sample vial depth, Injection depth
Plunger sampling and injection speed,
Solvent plug volume,
Internal standard addition,
Air plug volume


Sample/Internal standard rinse
Syringe strokes before injection
Pre- and post- injection delay
Viscosity delay
Sample volume

System robustness and internal design simplicity are combined to ideally perform routine and research laboratory workload providing fast, reliable, and accurate analysis.

DANI

The demand of increased productivity is more and more frequent in modern laboratories along with high reproducibility and repeatability requirements.

Fatty Acid Methyl Esters




The Master AS – Dual Injection allows two samples to be simultaneously injected into two different channels, therefore performing two analyses at the same time. The determination of FAMES demonstrates the Master AS Dual Injection performance. The run time of this analysis is about 55 minutes, therefore a laboratory can process around 25 samples a day. Thanks to Dual Injection Mode, 50 analyses have been run with excellent results in terms of area and retention time repeatability

#	FAME	RT Det A	RT Det B	RT St.Dev. Det A	RT St.Dev. Det B	Area % St.Dev. Det A	Area % St.Dev. Det B
1	Methyl Butyrate	4.717	4.780	0.0062	0.0032	3.11%	3.16%
2	Methyl Hexanoate	8.898	8.961	0.0035	0.0051	3.57%	2.87%
3	Methyl Octanoate	14.021	14.084	0.0044	0.0022	3.69%	2.99%
4	Methyl Decanoate	19.029	19.092	0.0051	0.0063	3.95%	3.53%
5	Methyl Undecanoate	21.376	21.439	0.0030	0.0071	3.96%	3.71%
6	Methyl Laurate	23.641	23.704	0.0028	0.0048	3.99%	3.61%
7	Methyl Tridecanoate	25.791	25.854	0.0045	0.0034	3.86%	3.52%
8	Methyl Tetradecanoate	27.876	27.939	0.0063	0.0035	3.85%	3.38%
9	Myristoleic Acid Methyl Ester	28.591	28.654	0.0074	0.0027	3.85%	3.38%
10	Methyl Pentadecanoate	29.849	29.912	0.0021	0.0060	3.81%	3.25%
11	cis-10-Pentadecenoic Acid Methyl Ester	30.554	30.619	0.0036	0.0055	3.81%	3.26%
12	Methyl Palmitate	31.777	31.842	0.0051	0.0036	3.73%	3.12%
13	Methyl Palmitoleate	32.229	32.294	0.0033	0.0029	3.57%	3.22%
14	Methyl Heptadecanoate	33.578	33.643	0.0076	0.0061	3.16%	2.95%
15	cis-10-Heptadecenoic Acid Methyl Ester	34.038	34.103	0.0064	0.0047	3.37%	2.94%
16	Methyl Octadecanoate	35.357	35.422	0.0033	0.0028	3.16%	3.18%
17	trans-9-Elaidic Acid Methyl Ester	35.702	35.767	0.0058	0.0066	3.37%	3.07%
18	cis-9-Oleic Acid Methyl Ester	35.702	35.767	0.0058	0.0066	3.37%	3.07%
19	Linoleic Acid Methyl Ester	36.454	36.519	0.0050	0.0029	3.00%	3.18%
20	Methyl Linoleate	36.517	36.582	0.0022	0.0055	3.34%	3.27%
21	Methyl Arachidate	36.962	37.027	0.0039	0.0079	3.19%	3.35%
22	gamma-Linolenic Acid Methyl Ester	37.523	37.596	0.0073	0.0082	3.23%	3.23%
23	Methyl cis-11-Eicosenoate	38.683	38.756	0.0071	0.0068	3.43%	3.32%
24	Methyl Linolenate	38.996	39.069	0.0057	0.0052	3.29%	3.17%
25	Methyl Heneicosanoate	39.741	39.814	0.0045	0.0049	3.31%	3.37%
26	cis-11,14-Eicosadienoic Acid Methyl Ester	40.167	40.240	0.0039	0.0067	3.83%	4.04%
27	Methyl Docosanoate	40.249	40.322	0.0036	0.0044	2.89%	3.17%
28	cis-8,11,14-Eicosatrienoic Acid Methyl Ester	40.540	40.613	0.0022	0.0026	3.35%	3.19%
29	Methyl Erucaate	40.788	40.861	0.0026	0.0020	3.34%	3.33%
30	cis-11,14,17-Eicosatrienoic Acid Methyl Ester	41.691	41.764	0.0034	0.0019	3.32%	3.05%
31	Methyl Tricosanoate	42.003	42.076	0.0071	0.0028	3.62%	3.64%
32	Methyl cis-5,8,11,14-Eicosatetraenoic	42.417	42.495	0.0069	0.0044	3.50%	2.92%
33	cis-13-16-Docosadienoic Acid Methyl Ester	43.423	43.501	0.0062	0.0055	3.64%	3.25%
34	Methyl Lignocerate	44.090	44.168	0.0043	0.0068	3.68%	3.59%
35	Methyl cis-5,8,11,14,17-Eicosapentaenoate	46.717	46.795	0.0044	0.0088	3.90%	3.49%
36	Methyl Nervonate	47.221	47.299	0.0038	0.0050	3.71%	3.80%
37	All cis-4-7-10-13-16-19-Docosahexaenoate	47.381	47.459	0.0067	0.0062	3.44%	3.92%

Nitro/phosphorous and Organochloride Pesticides

Chlorinated and nitrogen/phosphorus compounds are determined by two selective detectors (ECD and NPD).

These two different channels/detectors allow different classes of compounds to be simultaneously detected without affecting sensitivity like in systems where flow is split in two columns.



ECD [in red] and NPD [in blue] chromatogram

Nitro/Phosphorus Pesticides
Tionazine
Phorate
Diazinone
CH3-Parathion
Parathion
Quinalphos
Metidathion
Phosdrin
Phonophos
Pirimiphos

Organochloride Pesticides
Aldrin
alpha-BHC
beta-BHC
Lindane
delta-BHC
delta-BHC
gamma-Chlordane
4,4'-DDE
1,1-Dichloro-2,2-bis(4-chlorophenyl)ethane
4,4'-DDT
Decachlorobiphenyl
Dieldrin
alpha-Endosulfan
beta-Endosulfan
beta-Endosulfan
Endrin
Endrin Aldehyde
Endrin Ketone
Heptachlor
Heptachlor Exo-Epoide
Methoxychlor
2,4,5,6-Tetrachloro-m-xylene



AUTHORIZED DEALER

Disclaimer

The contents of this brochure are for reference and illustrative purposes only. Information, descriptions and specifications in this publication are subjected to change without notice.
DANI Analitica srl assumes no responsibility and will not be liable for any errors or omissions contained herein or for incidental, consequential damages or losses in connection with the furnishing, performance or use of this material.

For a complete listing of our products visit www.dani-analitica.com